

Lojer Field Operating Table

Tested by civilian and military field hospitals all around the world

Lojer Field Operating Table has been developed in close cooperation with the Red Cross and military medical experts, and it has been used in civilian and military field hospitals all around the world. The table also has Nato-coding.

Field Operating Table is manufactured from high quality raw materials. The metal frame is made of epoxy-coated steel and table top is made of waterproof plywood board. Table is fully mechanical; all the functions are manually operated, nevertheless providing sufficient number of features to execute efficient field operations. The table also features 10 x 25 mm siderails for accessories. Installation and operation are extremely easy and when the table is not in use it can be folded to save space.

The entire table is designed to fit into a carrying box measuring 78 x 38 x 41 cm (= 0,12 m³), and the table weight approx. 65 kg. The aluminium carrying box is built for rough conditions and is dimensioned according to the Red Cross standards. The mattress is 50 mm thick, size 190 x 50 cm (in four sections) and all the other accessories can be packed in a similar carrying box.

Ready for field operations in just few minutes.



Lojer Field Operating Table

LOJER® For easy care

Field Operating Table	
Carrying box - aluminium, weight 9 kg	78 x 38 x 41 cm
Overall length (with head plate)	2020 mm
Width	500 mm
Safe Working Load (SWL)	150 kg
Table weight	65 kg
Table height	680..900 mm (50 mm intervals)
Base Stand, aluminium 147 x 124 cm (3 x 49 cm x124 cm) with fixing levers	8 kg
Trendelenburg / Reverse Trendelenburg	20°/ 32°
Gas spring assisted back section (up/down)	30°/ 72°
Gas spring assisted leg section (up/down)	26°/ 83°
Head section (up/down), detachable	29°/ 87°
Warranty	2 years
Country of Origin	Finland

- Legs have a holes for fixing the table to the floor, ground or specially made aluminium base stand.
- Entire foot section can easily be removed and gynaecological section can be installed
- List of accessories:



1. IV-pole
2. Rapid rod clamp
3. Anaesthesia frame
4. Side supports
5. Gynaecological kit featuring fasteners, stainless steel bowl and pair of leg supports
6. Arm board

