

LOJERIN (JA MERIVAARAN) SÄNGYT

Sairaalasängyt

Malli:

ScanAfia X S/ TK/ ICU

leveys 80cm	78x205xmax.13cm
leveys 90cm	88x205xmax.13cm

Sallitut patjakoot erityisehdoin:

78x205xmax.21cm (korotetuin kaitein 70355KP)
88x205xmax.21cm (korotetuin kaitein 70355KP)

Patjavaihtoehdot:

FlexCare/ FlexCare Rescue
HighCare/ HighCare Rescue

ScanAfia X ICU W

leveys 80cm	78x205xmax.13cm
leveys 90cm	88x205xmax.13cm

78x205xmax.21cm (korotetuin kaitein 71034KP)
88x205xmax.21cm (korotetuin kaitein 71034KP)

FlexCare/ FlexCare Rescue
HighCare/ HighCare Rescue

Carena + W

leveys 80cm	75x200x12/13
leveys 85cm	80x200x12/13
leveys 90cm	85x200x12/13

FlexCare/ FlexCare Rescue
HighCare/ HighCare Rescue

Futura Plus

leveys 81,2 cm	80x200x12/13
----------------	--------------

FlexCare/ FlexCare Rescue
HighCare/ HighCare Rescue

Hoitosängyt

Malli:

Modux

leveys 80cm	78x205x(9)12-13cm
leveys 90cm	88x205x(9)12-13cm

78x205xmin.9-max.21cm (kaidekoroke 71025KP, 71026KP)
88x205xmin.9-max.21cm (kaidekoroke 71025KP, 71026KP)

FlexCare/ FlexCare Rescue
HighCare/ HighCare Rescue

ScanAfia X HS

leveys 80cm	78x205xmin.9cm-max.13cm
leveys 90cm	88x205xmin.9cm-max.13cm
leveys 80cm jatkettuna	78x215xmin.9cm-max.13cm
leveys 90cm jatkettuna	88x215xmin.9cm-max.13cm

78x205xmin.9cm-max.21cm (kaidekoroke 71025KP, 71026KP)
88x205xmin.9cm-max.21cm (kaidekoroke 71025KP, 71026KP)
Kaidekorotus ei tähän malliin mahdollinen
Kaidekorotus ei tähän malliin mahdollinen

FlexCare/ FlexCare Rescue
HighCare/ HighCare Rescue

ScanAfia Pro HS 280, 290, 480, 490, 480 215, 490 215

leveys 80cm	78x205xmin.9cm-max.13cm
leveys 90cm	88x205xmin.9cm-max.13cm
leveys 80cm jatkettuna	78x215xmin.9cm-max.13cm
leveys 90cm jatkettuna	88x215xmin.9cm-max.13cm

78x205xmin.9cm-max.21cm (kaidekoroke 71025KP, 71026KP)
88x205xmin.9cm-max.21cm (kaidekoroke 71025KP, 71026KP)
Kaidekorotus ei tähän malliin mahdollinen
Kaidekorotus ei tähän malliin mahdollinen

FlexCare/ FlexCare Rescue
HighCare/ HighCare Rescue

ScanAfia X XL

leveys 120cm	118cmx205xmin.9cm-max.13cm
--------------	----------------------------

118cmx205xmin.9cm-max.21cm (kaidekoroke 71025KP, 71026KP)

HighCare/ HighCare Rescue
MaxCare/ MaxCare Rescue