

ScanAfia X ICU L hospital bed

Also available in a machine-washable version



Main features:



1. Very high load capacity: SWL 340 kg
2. Large twin castors (150 mm) and lightweight design make the bed effortless to move around
3. Impressive 20° anti-Trendelenburg and Trendelenburg adjustment
4. Integrated mattress platform extension

The ScanAfia X ICU L hospital bed is highly suitable for a variety of hospital wards. The bed's exceptional height adjustment range improves ergonomics for healthcare professionals and supports patient self-sufficiency. In addition to a handset, the nurse control panel (supplied as standard) allows for securing the bed's functions so that the patient cannot use them. It also enables specialised electrical functions, including the Trendelenburg, CPR and cardiac positions. The sliding mechanism in the back section prevents the patient from slipping toward the foot-end of the bed and improves comfort. Quick-release locks at the head and foot of the bed allow for the rapid detachment of end frames in emergencies. Despite the substantial 340 kg safe working load (SWL), the standard 150 mm twin castors make ScanAfia X ICU L easy and effortless to move around.

The bed is designed for long-term use in hospital environment, and it is available with a width of 80 cm or 90 cm. A washable version of this hospital bed model is also available, complete with a wash cycle counter.

Technical data and features	
Total width (with siderails)	89.5 cm / 99.5 cm
Patient surface width	80 cm / 90 cm
Total length	217 cm (extended 245 cm)
Patient surface length	205 cm (extended 225 cm)
Height adjustment	34–84 cm
Recommended mattress size	78 / 88 x 205 x 13 cm
Total weight (without rails)	115 kg / 80 cm 117 kg / 90 cm
Castors	Ø150 mm central locking twin castors with directional wheel
Safe Working Load (SWL)	340 kg
Electrical system	24V
Back section adjustment	70°
Thigh section adjustment	30°
Calf section adjustment	18°
Trendelenburg / anti-Trendelenburg	20° / 20°
Removable bed ends	White, blue, yellow and light green
Side rails	Different options available
Other features	Battery backup hand control nurse control panel back section slide electrical CPR-function manual emergency lowering of back section fixing hooks for supporting belts roller bumpers
CE marked	Yes (The device is a CE marked medical device according to MDR 2017/745 class 1)
Warranty	2 years
Country of Manufacture	Finland (Lojer Oy)

Factory installed options:

- 5th castor for easier manoeuvring
- Various castor sizes available
- Wash cycle counter
- Led or lithium-ion battery
- Chrome-plated side rails and end frames
- Raised rails

Other accessories:

- Handset holder
- Grab handle
- Headboard
- IV stand
- Lifting pole
- Overbed dining table
- Urine bag holder
- Monitor and writing stand
- Oxygen cylinder holder
- Accessory mounting rail

